Sodium Iodide Symporter (NIS) Expression in Cytologically Indeterminate and Malignant Thyroid Nodules

Prasana Santhanam MBBS, MD,¹ Mohammed Alshalalfa, PhD,² Yangyang Hao, PhD,² Joshua P. Klopper, MD,² Richard T. Kloos, MD,² Paul W. Ladenson, MD¹ 1. The Johns Hopkins University School of Medicine, Baltimore, Maryland, USA. 2. Veracyte, Inc, South San Francisco, CA, USA.

INTRODUCTION

- Well differentiated thyroid cancers vary significantly in their ability to concentrate radioiodine and, consequently, be accurately imaged and effectively treated with ¹³¹iodine.
- The sodium iodide symporter (NIS) is a transmembrane glycoprotein (codified by the SLC5A5 gene) expressed predominantly in the basolateral plasma membrane of thyroid follicular cells.
- NIS mediates the active transport of iodine from the bloodstream into thyroid follicular cells, playing an essential role in the treatment and monitoring of differentiated thyroid carcinomas when radioiodine therapy is indicated.
- The objective of this study was to evaluate the preoperative mRNA expression of NIS in FNA samples of cytologically indeterminate (Cyto I – Bethesda (B) III or IV) and suspicious for malignancy or malignant thyroid (B V and VI) nodules sent for molecular testing and associate its expression with genomic alterations.

METHODS

- NIS mRNA expression was analyzed in 47,695 thyroid nodules sent for Afirma Genomic Sequencing Classifier (GSC) molecular testing.
- Differential NIS expression was explored across Afirma benign (GSC-B) and suspicious (GSC-S) categories, Bethesda cytology categories (B III-VI), and relative to various molecular alterations identified by the Afirma GSC and Xpression Atlas.

RESULTS

- Relative to B III/IV GSC-B nodules, B III/IV GSC-S and B V/VI nodules had lower NIS expression with a mean fold reduction (FC) 0.86 and 0.80 respectively (p < 0.01 for both) (Figure 1).
- Among B III/IV GSC-S nodules (Figure 2a) and B V/VI nodules (Figure 2b), NIS expression was different based upon the presence of certain molecular alterations detected by Afirma XA. NIS expression was lowest in nodules with BRAFV600E variant expression and highest in nodules with TSHR variant expression.
- Nodules with a high Hürthle cell index (a component of the Afirma GSC core ensemble classifier) and high-risk molecular features of Hürthle cell carcinoma had lower NIS expression (Figures 3a and 3b).



Within GSC-S group, NIS was lowest in samples with BRAFV600E and highest in samples with Thyroid Stimulating Hormone Receptor (TSHR). Compared to GSC-B samples, nodules with BRAFV600E or RAS variants had lower expression (FC:0.80 and 0.84 respectively, p<2e⁻¹⁶ for both), samples with TSHR variants had higher expression (FC:1.49, p<2e⁻¹⁶), and SPOP variants and *RET/ALK* fusion were not different (FC:0.96 and 0.93 respectively).

When evaluating NIS expression with relation to BRAFV600 in the GSC-S group, NIS expression was highest in samples with Thyroid Stimulating Hormone Receptor (TSHR) (FC:1.87 relative to BRAFV600E samples, p<2e⁻¹⁶) variants, SPOP (FC:1.2, p=5e⁻¹⁰), and *RET/ALK* fusion (FC:1.17, p<2e⁻¹⁶). (Data not shown on graph.)



Within Bethesda V/VI nodules, NIS was lowest in samples with BRAFV600E and highest in samples with Thyroid Stimulating Hormone Receptor (*TSHR*). Relative to GSC-B samples, Bethesda V/VI nodules harboring *BRAF*V600E had lower NIS expression (FC:0.73, p<2e⁻¹⁶), whereas samples with TSHR variants had higher expression (FC:1.34, p=0.003) and SPOP samples had no difference (FC:1, p=0.9).

FIGURE 3a



FIGURE 3b



CONCLUSION

- NIS expression in indeterminate and malignant thyroid nodules varies across Bethesda cytology groups, Afirma GSC classification, and certain molecular alterations.
- NIS expression is necessary and may not be sufficient for thyroid tumor response to radioiodine.¹
- To address the clinically relevant predictive significance of NIS mRNA expression levels in preoperative biopsy specimens, future studies should evaluate tumor iodine avidity and the response to radioiodine therapy.

REFERENCE





In BIII/IV GSC-S nodules, there is an inverse correlation between the Afirma Hürthle cell index and NIS expression across all tested nodules. Note that the red rectangle highlights that almost all samples with a Hürthle high index score have low NIS mRNA expression.

> BIII/IV GSC-S samples with high Hürthle index, high neoplasm index, and high loss of heterozygosity (LOH) are more likely to be Hürthle cell carcinoma and these samples have the lowest NIS expression.